

# Catalog Number

68340-25P-4

**HPLC Column, NP Amino, 3 $\mu$ m, 4.0mm ID x 250mm Length, 100A.**  
Cogent Brand. 1 EA. CUSTOM ORDER.



## Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 3 $\mu$ m Particle, this column is more efficient than a 5 $\mu$ m column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure.

**Price****\$1,036.93**

## Additional Info

Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,036.93**