

Catalog Number

68340-25P-4-3M

HPLC Column, NP Amino, Wide Pore, 3 μ m, 4.0mm ID x 250mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent NP Amino™ Wide Pore Stationary Phase has been a value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group, to the highest purity silica available. This makes this column ideal as an Alternate in USP Methods that use the L8 Classification, for large molecules.

With a 3 μ m Particle and 300A Pore Size, this column performs well for polar molecules with a negative charge such as complex Sugars and larger compounds with Carboxylic Acid Groups using Normal Phase or *HILIC* modes.

Price**\$1,036.93**

Additional Info

These columns are packed and shipped in Normal Phase Solvents.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,036.93