

Catalog Number

68504-15P-4-3M

HPLC Column, RP C4 Butyl, Wide Pore, 5 μ m, 4.0mm ID x 150mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available.

With a 5 μ m Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Used for proteins, peptides, enzymes and other large, hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for bioanalysis.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$848.36

Printed On 2024/08/06

Additional Info

order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$848.36

Printed On 2024/08/06