

# Catalog Number

68504-P30-100

**Prep HPLC Column, RP C4, 5um, 100A. 30mm ID x 100mm Length.**  
Cogent Brand. 1 EA.



## Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 carbon (*butyl*) ligand to the highest purity silica available.

The 5 µm Particle works well with prep flow rates and produces lower back pressures than smaller particles. The large diameter and length of this column provides high loading capacity which is great for purification of hydrophobic or neutral compounds making this it ideal for Preparative HPLC.

**Price****\$6,843.03**

## Additional Info

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$6,843.03**