

## Catalog Number

68504-P30-2503M

**Prep HPLC Column, RP C4, Wide Pore, 5 $\mu$ m, 300A. 30mm ID x 250mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available. This makes this column ideal for large molecules.

The 5  $\mu$ m Particle works well with prep flow rates and produces lower back pressures than smaller particles. The large diameter and length of this column provides high loading capacity which is great for purification of proteins, peptides, enzymes and other large, hydrophobic or neutral compounds in Reversed Phase Preparative HPLC.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$11,612.43**

Printed On 2024/05/17