

## Catalog Number

68504-P50-50-3M

**Process HPLC Column, RP C4, Wide Pore, 5µm, 300A. 50mm ID x 50mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available. This makes this column ideal for large molecules.

With a 5µm Particle and 300A Pore Size, this column works well with high flow rates and acceptable back pressures. Used in process for purification of proteins, peptides, enzymes and other large, hydrophobic or neutral compounds in Reversed Phase preparative HPLC, this Butyl Phase is excellent for bio molecules.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$10,914.16**

Printed On 2024/08/06

## Additional Info

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$10,914.16**

Printed On 2024/08/06