

# Catalog Number

68504-PG1-3

**Guard Column Cartridges, RP C4 Wide Pore, 5µm, 300A. 2.1mm ID x 10mm Length. Cogent Brand. 5/PK.**



## Additional Info

These Guard Column Cartridges are packed with 5µm, Cogent RP C4™ 300A, Wide Pore Stationary Phase. They are double end capped (*fully*), with traditional bonding of a 4 carbon ligand (Butyl) to the highest purity silica available.

Designed to protect your 2.1mm ID, RP C4, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder,*

**Price****\$458.03**

## Additional Info

click [HERE](#) for more information about the Holder.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$458.03**