

## Catalog Number

68504-SP150-3M

**Semi Prep HPLC Column, RP C4 Wide Pore, 5um, 300A. 10mm ID x 150mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C4™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*Butyl*) carbon ligand to the highest purity silica available.

With a 5µm Particle and 300A Pore Size, this column works well with high flow rates and acceptable back pressures. The diameter and length of this column provides increased loading capacity compared with analytical columns, which is why it is chosen for purification of proteins, peptides, enzymes and other biological compounds making it ideal for Semi Preparative HPLC.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,277.36**

Printed On 2024/05/18