

## Catalog Number

68508-25P-3M

**HPLC Column, RP C8 Octyl, Wide Pore, 5 $\mu$ m, 4.6mm ID x 250mm Length, 300A. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C8™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available.

With a 5  $\mu$ m Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for large hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is also excellent for long chain fatty acids, some proteins and peptides.

*This column is often chosen for use macro-molecules such as proteins and polysaccharides in HPLC using UV or ELSD detection.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$921.70**

Printed On 2024/05/02