

Catalog Number

68508-P21-2503M

Prep HPLC Column, RP C8 Octyl, Wide Pore, 5um, 300A. 21.2mm ID x 250mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent RP C8™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories and pilot plants for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available.

The 5 µm Particle works well with prep flow rates and produces lower back pressures than smaller particles. The large diameter and extra length of this column provides high loading capacity which is great for purification of large compounds making it ideal for many proteins & peptides and other compounds.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$6,230.53

Printed On 2024/04/28