

Catalog Number

68512-25P-4

HPLC Column, RP C12, 5um, 4.0mm ID x 250mm Length, 100A.
Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent RP C12™ Stationary Phase has been a premium value and used by HPLC laboratories for several years. It is double end capped (*fully*), with bonding of a rigid, Adamantine shaped carbon silane (diamond like) to the highest purity silica available, making this column excellent for complex mixtures with rapid equilibration between runs.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Compatible with any Reversed Phase mobile phase including 100% water, you will not experience any loss of retention run

Price**\$1,070.95**

Additional Info

to run and therefore this column is often used for high throughput screening (HTS) methods. Normally, you can expect less retention than a C18 column but more than a C8 column for the same hydrophobic compounds.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,070.95