

Catalog Number

68512-P21-250

Prep HPLC Column, RP C12, 5um, 100A. 21.2mm ID x 250mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP C12™ Stationary Phase has been a premium value and used by HPLC laboratories for several years. It is double end capped (*fully*), with bonding of a rigid, Adamantine shaped carbon silane (diamond like) to the highest purity silica available, making this column excellent for complex mixtures with rapid equilibration between runs.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Compatible with any Reversed Phase mobile phase including 100% water, you will not experience any loss of retention run to run and therefore this column is often used for preparative HPLC methods. Normally, you can expect less retention than a C18 column but more than a C8 column for the same hydrophobic compounds.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$6,950.26

Printed On 2024/07/27