

Catalog Number

68512-P21-PG1

Prep Guard Column Cartridge, RP C12, 5um, 100A. 21.2mm ID x 10mm Length. Cogent Brand. 1 EA.



Additional Info

This Guard Column Cartridge is packed with 5µm, Cogent RP C12™ Stationary Phase. It is double end capped (*fully*), with bonding of a rigid, Adamantine shaped carbon silane (diamond like) to the highest purity silica available. It is designed to protect your RP C12 Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

Use this Guard Column Cartridge with a 21.2mm ID Prep Column Holder. Click [HERE](#) for more information.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$447.95

Printed On 2024/05/08