

# Catalog Number

68512-PG2

**Guard Column Cartridges, RP C12, 5 $\mu$ m, 100A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.**



## Additional Info

These Guard Column Cartridges are packed with 5 $\mu$ m, Cogent RP C12™ Stationary Phase. They are double end capped (*fully*), with bonding of a rigid, Adamantine shaped carbon silane (diamond like) to the highest purity silica available.

Designed to protect your 3.0, 4.0 or 4.6mm ID, RP C12, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder,*

**Price****\$458.03**

## Additional Info

click [HERE](#) for more information about the Holder.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$458.03**