

## Catalog Number

68518-05P-3M

**HPLC Column, RP C18, Wide Pore, 5um, 4.6mm ID x 50mm Length, 300A. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C18™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available.

With a 5 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for large hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is also excellent for long chain fatty acids, some proteins and peptides.

*This size column provides good resolution and short run times and is Very Popular Column used in High Through-put Screening.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$733.16**

Printed On 2024/07/27