

Catalog Number

68518-15P-3M

HPLC Column, RP C18, Wide Pore, 5µm, 4.6mm ID x 150mm Length, 300A. Cogent Brand. 1 EA.



Additional Info

This Cogent RP C18™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available.

With a 5 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures. Used for large hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is also excellent for long chain fatty acids, some proteins and peptides.

This column is often chosen for use macro-molecules such as proteins and polysaccharides in HPLC using UV detection or in LCMS with a flow splitter.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$848.36

Printed On 2024/04/25