

Catalog Number

68518-P30-100

Prep HPLC Column, RP C18, 5um, 100A. 30mm ID x 100mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP C18™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available.

The 5 µm Particle works well with prep flow rates and produces lower back pressures than smaller particles. The extra large diameter and long length of this column provides very high loading capacity which is great for purification of hydrophobic or neutral compounds making this it ideal for Preparative HPLC.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$6,692.45

Printed On 2024/08/06

Additional Info

ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$6,692.45

Printed On 2024/08/06