

## Catalog Number

68518-SP-PG1-3

**Semi Prep Guard Column Cartridge, RP C18 Wide Pore, 5 $\mu$ m, 300A.**  
10mm ID x 10mm Length. Cogent Brand. 5/PK.



### Additional Info

These Guard Column Cartridges are packed with 5 $\mu$ m, Cogent RP C18™ Wide Pore Stationary Phase. They are double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available. They are designed to protect your RP C18 Wide Pore Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in Stainless Steel Tubes, they have Stainless Steel Frits.

*Use these Guard Column Cartridges with a 10mm ID Prep Column Holder. Click [HERE](#) for more information.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$895.91**

Printed On 2024/05/02