

# Catalog Number

68518-SP150

**Semi Prep HPLC Column, RP C18, 5 $\mu$ m, 100A. 10mm ID x 150mm Length. Cogent Brand. 1 EA.**



## Additional Info

This Cogent RP C18™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available.

The 5 $\mu$ m Particle works well with high flow rates and produces lower back pressures than smaller particles. The diameter and length of this column provides increased loading capacity compared with analytical columns, which is why it is chosen for purification of hydrophobic or neutral compounds, making it ideal for Semi Preparative HPLC.

## Additional Info

**Notes** - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)