

Catalog Number

68522-10P-4

HPLC Column, RP PFP, 5µm, 4.0mm ID x 100mm Length, 100A.
Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines; the highest purity silica available is used making this column ideal as an Alternate in USP Methods that require the L43 Classification.

With a 5µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Often chosen for samples with halogenated or polar functional groups.

Packed and shipped in Reversed Phase Solvents.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$785.53

Printed On 2024/04/28