

## Catalog Number

68522-25P-4

**HPLC Column, RP PFP, 5µm, 4.0mm ID x 250mm Length, 100A.**  
Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines; the highest purity silica available is used making this column ideal as an Alternate in USP Methods that require the L43 Classification.

With a 5µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Often chosen for samples with halogenated or polar functional groups.

*Packed and shipped in Reversed Phase Solvents.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$921.70**

Printed On 2024/05/08