

Catalog Number

68522-SP-PG1

Semi Prep Guard Column Cartridge, RP PFP, 5 μ m, 100A. 10mm ID x 10mm Length. Cogent Brand. 5/PK.



Additional Info

These Guard Column Cartridges are packed with 5 μ m, Cogent RP PFP™ Pentafluorophenyl Stationary Phase. It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

They are designed to protect your RP PFP™ Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, they have Stainless Steel Frits.

Use these Guard Column Cartridges with a 10mm ID Prep Column Holder. Click [HERE](#) for more information.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$895.91

Printed On 2024/05/05