

## Catalog Number

68524-P21-50-3M

**Prep HPLC Column, RP Phenyl, Wide Pore, 5 $\mu$ m, 300A. 21.2mm ID x 50mm Length. Cogent Brand. 1 EA.**



### Additional Info

The Cogent RP Phenyl™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories and pilot plants for many years. It is double end capped (fully), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 5 $\mu$ m particle with 300A pore size works well with higher flow rates.

The large diameter and short length of this column provides good loading capacity which is great for purification of larger compounds such as isomers (positional),

## Additional Info

Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Preparative HPLC.

**Notes** - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)