

Catalog Number

68524-P30-50

Prep HPLC Column, RP Phenyl, 5um, 100A. 30mm ID x 50mm Length.
Cogent Brand. 1 EA.



Additional Info

The Cogent RP Phenyl™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 5µm particle works well with higher flow rates.

The extra large diameter and short length of this column provides good loading capacity which is great for purification of isomers (positional), Tocopherols,

Additional Info

Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Preparative HPLC.

Notes - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com