

## Catalog Number

68524-P50-150

**Process HPLC Column, RP Phenyl, 5 $\mu$ m, 100A. 50mm ID x 150mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 5 $\mu$ m particle works well with higher flow rates.

The extra large diameter and length of this column provides high loading capacity which is great for purification of isomers (positional), Tocopherols, Flavonoids,

## Additional Info

Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Process HPLC.

**Notes** - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)