

Catalog Number

68524-P50-2503M

Process HPLC Column, RP Phenyl, Wide Pore, 5µm, 300Å. 50mm ID x 250mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 5µm particle with 300Å pore size works well with higher flow rates.

The extra large diameter and very long length of this column provides very high loading capacity which is great for purification of larger compounds such as isomers

Additional Info

(positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Process HPLC.

Notes - *A specialized HPLC System is required for columns this size.* This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com