

# Catalog Number

68524-SP-PG1-3

**Semi Prep Guard Column Cartridge, RP Phenyl, Wide Pore, 5µm, 300A. 10mm ID x 10mm Length. Cogent Brand. 5/PK.**



## Additional Info

These Guard Column Cartridges are packed with 5µm, Cogent RP Phenyl™ Wide Pore, Stationary Phase. They are double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available. They are designed to protect your 5µm RP Phenyl Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, they have Stainless Steel Frits.

*Use these Guard Column Cartridges with a 10mm ID Prep Column Holder.*

*Click [HERE](#) for more information.*

**Price****\$916.07**

## Additional Info

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$916.07**