

## Catalog Number

68534-25P

**HPLC Column, RP Cyano**, 5 $\mu$ m, 4.6mm ID x 250mm Length, 100A.

Cogent Brand. 1 EA.

**MOST POPULAR.**



### Additional Info

This Cogent RP Cyano™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond, and it is bonded to the highest purity silica available making this Column ideal as an Alternate in USP Methods that use the L10 Classification.

With a 5  $\mu$ m Particle, this Column works well with higher flow rates and lower back pressures compared with smaller Particles. It is often used for polar or neutral compounds in Reversed Phase HPLC and you should expect shorter retention times and faster equilibration (*between runs*) than with a non-modified, Silica Column.

Packed and shipped in Reversed Phase Solvents; click [HERE](#) for Cyano Columns packed in Normal Phase (non polar) Solvents. *This is one of our Most Popular RP Cyano Columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$921.70**

Printed On 2024/04/28