

# Catalog Number

68534-PG2-3

**Guard Column Cartridges, RP Cyano Wide Pore, 5um, 300A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.**



## Additional Info

These Guard Column Cartridges are packed with 5µm, Cogent RP Cyano™ Wide Pore Stationary Phase. Not end capped, they are manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond and are bonded to the highest purity silica available.

Designed to protect your 3.0, 4.0 & 4.6mm ID, RP Cyano, Analytical Columns from particulates and contaminants that will adsorb to it. They are packed with Reversed Phase Solvent in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

**Price****\$458.03**

## Additional Info

*These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.*

### **MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### **Price**

**\$458.03**