

Catalog Number

68539-15P

HPLC Column, RP Phenyl Hexyl, 5um, 4.6mm ID x 150mm Length, 100A. Cogent Brand. 1 EA.

MOST POPULAR



Additional Info

This Cogent RP Phenyl Hexyl™ HPLC Column has been a premium value and used by HPLC laboratories for years. It is double end capped (*fully*), with traditional bonding of a 6 carbon silane with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification with Phenyl Hexyl stationary phase or for new methods.

Used to Separate aromatic compounds that have some hydrophobicity, these columns may offer different Selectivity from Phenyl Columns due to the extended C6 linker in the ligand. The pi-pi interactions of the Phenyl Group are heightened due to the distance from the Silica Surface (away from Silanol activity) and allows for excellent Retention and Separation of many polar compounds with Benzene rings. These include some

Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds. To optimize the pi interaction, using Methanol instead of Acetonitrile is suggested.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/04/28