

Catalog Number

68539-15P-2

HPLC Column, RP Phenyl Hexyl, 5um, 2.1mm ID x 150mm Length, 100A. Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl Hexyl™ HPLC Column has been a premium value and used by HPLC laboratories for years. It is double end capped (*fully*), with traditional bonding of a 6 carbon linkage with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in some USP Methods that use the L11 Classification and a Phenyl Hexyl Stationary Phase.

Used to Separate aromatic compounds that have some hydrophobicity, these columns may offer different Selectivity from Phenyl Columns due to the extended C6 linker in the ligand. The pi-pi interactions of the Phenyl Group are heightened due to the distance from the Silica Surface (away from Silanol activity) and allows for excellent Retention and Separation of many polar

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$995.01

Printed On 2024/08/06

Additional Info

compounds with Benzene rings. These include some Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds. To optimize the pi interaction, using Methanol instead of Acetonitrile is suggested.

This column size is used with LCMS.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$995.01

Printed On 2024/08/06