

Catalog Number

68540-25P

HPLC Column, NP Amino, 5um, 4.6mm ID x 250mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

Elevate Your Research with Cogent NP Amino™ Stationary Phase! Trusted by Scientists | Ideal for USP L8 Methods | Versatile Applications

The Cogent NP Amino™ Stationary Phase has been a cornerstone in HPLC laboratories for years, offering exceptional value and performance. This column is not end-capped and features a 3-carbon ligand (propyl) with a terminal amino group bonded to world-class silica. These qualities make it an excellent alternative for any USP Method using the L8 Classification including the Grace Alltima™ Amino columns or for developing new methods.

Price**\$942.44**

Additional Info

Why Choose the Cogent NP Amino™?

- **High-Quality Silica:** Ensures reliable and reproducible results.
- **Versatile Applications:** Perfect for polar molecules with a negative charge, such as sugars and compounds with carboxylic acid groups, in Normal Phase or HILIC modes.
- **Optimal Performance:** Packed and shipped in Normal Phase solvents, this column is designed to deliver consistent, high-quality separations.

Key Benefits:

- **Robust and Reliable:** Known for its robustness, this column is a popular choice for both established USP methods and innovative new method development.
- **Scientifically Proven:** Trusted by scientists for its excellent performance and versatility in various applications.

Join the Experts

Many researchers rely on the Cogent NP Amino™ for its superior performance and reliability. Whether you're working with polar molecules or developing new methods, this column is designed to meet your scientific needs.

Ready to Enhance Your Lab's Data Quality?

Experience the difference with the **Cogent NP Amino™ Stationary Phase** and take your HPLC analyses to the next level!

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$942.44