

# Catalog Number

68540-P30-50-3M

**Prep HPLC Column, NP Amino, Wide Pore, 5um, 300A. 30mm ID x 50mm Length. Cogent Brand. 1 EA.**



## Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal for preparative HPLC.

Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

**Price****\$4,511.75**

## Additional Info

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$4,511.75**