

Catalog Number

68540-P30-PG1-3

Prep Guard Column Cartridge, NP Amino Wide Pore, 5 μ m, 300Å.
30mm ID x 10mm Length. Cogent Brand. 1 EA.



Additional Info

This Guard Column Cartridge is packed with 300Å, 5 μ m Cogent NP Amino™ Stationary Phase. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group, to the highest purity silica available. It is designed to protect your NP Amino Wide Pore, Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

Use this Guard Column Cartridge with a 30mm ID Prep Column Holder.
Click [HERE](#) for more information.

Additional Info

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com