

# Catalog Number

68540-SP250

**Semi Prep HPLC Column, NP Amino, 5um, 100A. 10mm ID x 250mm Length. Cogent Brand. 1 EA.**



## Additional Info

This Cogent NP Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group bonded to world class silica. This column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

Used for polar molecules with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or *HILIC* modes these columns are packed and shipped in Normal Phase Solvents.

**Price****\$2,674.84**

## Additional Info

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,674.84**