

Catalog Number

69020-P21-100

Prep HPLC Column, Phenyl Hydride, 4um, 100A. 21.2mm ID x 100mm Length, Preparative. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Phenyl Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for purification of compounds with aromatic rings that are difficult to retain or separate in Reversed Phase, including aromatic amines, small peptides and azides.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$7,489.80

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$7,489.80

Printed On 2024/08/06