

Catalog Number

69069-15P

HPLC Column, UDC Cholesterol, 4um, 4.6mm ID x 150mm Length, 100A. Cogent TYPE-C Brand. 1 EA.



Additional Info

This UDC Cholesterol™ HPLC is made with [Cogent TYPE-C™ Silica](#), and provides unique Selectivity due to different surface chemistry, compared to typical Reversed Phase Columns. This Silica Hydride Phase is modified with an 11 carbon chain with a cholesterol molecule attached at the terminus. This gives this column the ability to form “slots” and separate some compounds based on shape. Also, due to the direct silicon carbon bonds, this column has no need for end-capping. Use in Reversed Phase, Normal Phase and Aqueous Normal Phase (ANP).

Chosen for analyzing steroidal and other hard to retain compounds including some lipids, peptides, proteins & oligonucleotides.

- Faster Equilibration than other Reversed Phase or HILIC Columns
- Very Hydrophobic
- Robust Column Lifetimes

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,205.89

Printed On 2024/05/24