

## Catalog Number

69069-15P-2

**HPLC Column, UDC Cholesterol, 4um, 2.1mm ID x 150mm Length, 100A. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This UDC Cholesterol™ HPLC is made with [Cogent TYPE-C™ Silica](#), and provides unique Selectivity due to different surface chemistry, compared to typical Reversed Phase Columns. This Silica Hydride Phase is modified with an 11 carbon chain with a cholesterol molecule attached at the terminus. This gives this column the ability to form “slots” and separate some compounds based on shape. Also, due to the direct silicon carbon bonds, this column has no need for end-capping. Use in Reversed Phase, Normal Phase and Aqueous Normal Phase (ANP).

Chosen for analyzing steroidal and other hard to retain compounds including some lipids, peptides, proteins & oligonucleotides.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,205.89**

Printed On 2024/08/06

## Additional Info

- Faster Equilibration than other Reversed Phase or HILIC Columns
- Very Hydrophobic
- Robust Column Lifetimes
- Temperature and Solvent Changes Selectivity

*This column is often chosen for LCMS when increased resolution is desired.*

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$1,205.89**

Printed On 2024/08/06