

## Catalog Number

69220-10P-2

**HPLC Column, Phenyl Hydride 2.0**, 2.2 $\mu$ m, 2.1mm ID x 100mm  
Length, 120A, High Pressure. Cogent TYPE-C Brand. 1/EA



### Additional Info

This Cogent Phenyl Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for purification of compounds with aromatic rings that are difficult to retain or separate in Reversed Phase, including aromatic amines, small peptides and azides.

*This column is one of the more popular 2.0 column sizes and often chosen for LCMS when increased resolution is desired.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,295.33**

Printed On 2024/04/29