

## Catalog Number

## 69269-HG2

**Guard Column Starter Set, UDC Cholesterol 2.0.** Complete with 5 each 2.0mm ID x 10mm Long Cartridges with 2.2 $\mu$ m, 120A and 1 each Guard Column Holder. Cogent TYPE-C Brand. 1 PK.



### Additional Info

Everything you need to get started using Cogent UDC Cholesterol 2.0™ Guard Columns with your methods. This kit includes one "Universal" Stainless Steel [Guard Cartridge Column Holder](#) that can be used with any Cogent TYPE-C™ Guard Cartridge and five Cartridges, packed with Cogent UDC Cholesterol 2.0™ Stationary Phase are included.

This kit will normally be sufficient protection for the lifetime of an analytical column, based on average use. Each Cartridge in this Kit is shipped in individual, recloseable Delrin Cases to keep them from drying out or being damaged.

*Use this Starter Kit with 2.1mm ID analytical columns packed with Undecanoic Cholesterol 2.0 (UDC).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$1,478.80**

Printed On 2024/07/20