

Catalog Number

69269-PG3

Guard Column cartridge, UDC Cholesterol 2.0, replacement, 2.0 mm ID x 10 mm long cartridge packed with 2.2 μ m, 120 A phase. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, Delrin case to protect and keep it from drying out. Use with [guard column holder](#) sold separately.

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

Price**\$245.77**

Additional Info

Use this guard column cartridge with a 2.1 mm ID analytical column.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$245.77