

## Catalog Number

70000-10P-3

**HPLC Column, Diamond Hydride**, 4um, 3.0mm ID x 100mm Length,  
100A. Cogent TYPE-C Brand. 1 EA.

**Additional Info**

This Cogent Diamond Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites.

*This size column provides similar run times as a 4.6mm ID column, but uses much*

**Price****\$1,064.53**

## Additional Info

*less mobile phase, with lower flow rates. No special equipment is needed for this column to save solvent use.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,064.53**