

## Catalog Number

70000-15P-3

**HPLC Column, Diamond Hydride, 4um, 3.0mm ID x 150mm Length, 100A. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

**Unlock Superior Performance with Cogent Diamond Hydride HPLC Columns!**

Crafted with advanced Cogent [TYPE-C™ Silica](#), our Diamond Hydride HPLC Columns offer excellent selectivity, faster equilibration, and robust column lifetimes. Elevate your analytical capabilities with this innovative silica hydride-based stationary phase.

**Why Choose Cogent Diamond Hydride?**

**Price**

**\$1,233.02**

## Additional Info

- **Optimal for Polar Compounds:** Perfect for analyzing polar and non-polar peptides, amino acids, and other challenging biological compounds like metabolites.
- **Efficient and Economical:** Enjoy similar run times to a 4.6mm ID column while using much less mobile phase and lower flow rates. No special equipment needed, making it a cost-effective choice for solvent savings.

Custom Order Excellence: Each column is custom-manufactured to meet your specific needs. Please note that custom orders cannot be returned or canceled once placed.

**Ready to Transform Your Lab?** Order your Cogent Diamond Hydride HPLC Column today and experience unparalleled performance and reliability. Click [HERE](#) for more information on custom orders.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,233.02**