

## Catalog Number

70000-7.5P

**HPLC Column, Diamond Hydride, 4um, 4.6mm ID x 75mm Length, 100A. Cogent TYPE-C Brand. 1 EA. CUSTOM ORDER.**



### Additional Info

This Cogent Diamond Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites. When using in Aqueous Normal Phase, this column is Mass Spec Friendly.

This is a custom order, manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$941.74**

Printed On 2024/05/06