

## Catalog Number

70000-SP100

**Semi Prep HPLC Column, Diamond Hydride, 4um, 100A. 10mm ID x 100mm Length. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Diamond Hydride "Semi" Preparative, HPLC Column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for Purification or "Prep" of Peptides, Metabolites, Polar Compounds and other Analytes with amine groups when a large volume of sample is desired to be collected.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,708.86**

Printed On 2024/05/07