



Catalog Number

70200-15P-2

HPLC Column, Diamond Hydride 2.o, 2.2um, 2.1mm ID x 150mm Length, 120A. Cogent TYPE-C Brand. 1 EA.



Additional Info

Unlock Superior Performance with Cogent Diamond Hydride 2.o HPLC Columns!

Crafted with advanced Cogent [TYPE-C™ Silica](#), our Diamond Hydride 2.o HPLC Columns feature a smaller particle, silica hydride-based stationary phase that delivers increased efficiency, excellent selectivity, and faster equilibration. Achieve “Near UHPLC” performance on any HPLC instrument and elevate your analytical capabilities!

Price

\$1,601.51



Additional Info

Why Choose Cogent Diamond Hydride 2.0?

- **Superior Analysis of Polar Compounds:** Perfect for polar and non-polar peptides, amino acids, and other challenging biological compounds like metabolites.
- **Mass Spec Friendly:** Ideal for Aqueous Normal Phase applications, seamlessly integrating with your mass spectrometry workflows.
- **LCMS Ready:** Achieve increased resolution and precision in your LCMS analyses.

Ready to Transform Your Lab? Order your Cogent Diamond Hydride 2.0 HPLC Column today and experience unparalleled performance and reliability.

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,601.51