

Catalog Number

70200-15P-3

HPLC Column, Diamond Hydride 2.o, 2.2 μ m, 3.0mm ID x 150mm Length, 120A. Cogent TYPE-C Brand. 1 EA.



Additional Info

Unlock Superior Performance with Cogent Diamond Hydride 2.o HPLC Columns!

Crafted with advanced Cogent [TYPE-C™ Silica](#), our Diamond Hydride 2.o HPLC Columns feature a smaller particle, silica hydride-based stationary phase that delivers increased efficiency, excellent selectivity, and faster equilibration. Achieve “Near UHPLC” performance on any HPLC instrument and elevate your analytical capabilities!

Price

\$1,659.02

Additional Info

Why Choose Cogent Diamond Hydride 2.o?

- **Superior Analysis of Polar Compounds:** Perfect for polar and non-polar peptides, amino acids, and other challenging biological compounds like metabolites.
- **Efficient and Economical:** Enjoy similar run times to a 4.6mm ID column while using less mobile phase, lower flow rates, and experiencing less back pressure. No special equipment needed!

Custom Order Excellence: Each column is custom-manufactured to meet your specific needs. Please note that custom orders cannot be returned or canceled once placed.

Ready to Transform Your Lab? Order your Cogent Diamond Hydride 2.o HPLC Column today and experience unparalleled performance and reliability. Click [HERE](#) for more information on custom orders.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,659.02