

# Catalog Number

72030-PG2

**Guard Column Cartridges, C30, 10um, 200A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.**



## Additional Info

These Guard Column Cartridges are packed with 10µm, Cogent C30™ Stationary Phase. Not end capped, it is manufactured with traditional bonding of a ligand that has an average length of 30 carbons, to the highest purity silica available.

Designed to protect your 3.0, 4.0 or 4.6mm ID C30 Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder,*

**Price****\$458.03**

## Additional Info

click [HERE](#) for more information about the Holder.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$458.03**