

## Catalog Number

## 75020-P250

**Prep HPLC Column, HPS Phenyl, 5um, 21.2mm ID x 250mm Length, 120A. Cogent Brand. 1 EA.**



### Additional Info

This Cogent HPS Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 4 carbon ligand that has an attached *Phenyl* group and it is bonded to the highest purity, 120A silica available.

With a 5µm Particle, this Column works well with higher flow rates producing lower back pressures compared to columns with smaller particles making ideal for the purification of Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

## \$12,830.53

Printed On 2024/08/06

## Additional Info

*Compare this column's performance to Luna™ or Inertsil™ columns.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$12,830.53**

Printed On 2024/08/06