

## Catalog Number

## 75020-SP100

**Semi Prep HPLC Column, HPS Phenyl, 5um, 10mm ID x 100mm Length, 120A. Cogent Brand. 1 EA.**



### Additional Info

This Cogent HPS Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 4 carbon ligand that has an attached *Phenyl* group and it is bonded to the highest purity, 120A silica available.

With a 5µm Particle, this Column works well with higher flow rates producing lower back pressures compared to columns with smaller particles making ideal for the purification of Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC.

*Compare this column's performance to Luna™ or Inertsil™ columns.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

## \$2,986.04

Printed On 2024/05/05